

[index](#)**app** [/home/uwu/Github_Repos/University/M.O.S.I.S/M.O.S.I.S_Host_Software/App/app.py](#)

Launches app.py using python3 without calling interpreter explicitly.

Modules[os](#)**Classes**[builtins.object](#)[MediaEntryInternalRepresentation](#)[MediaMetadataInternalRepresentation](#)[enum.Enum\(builtins.object\)](#)[illuminationType](#)[shotType](#)sqlalchemy.orm.decl_api.Model([flask_sqlalchemy.model.Model](#))[MediaEntry](#)[MediaMetadata](#)class **MediaEntry**(sqlalchemy.orm.decl_api.Model)[MediaEntry](#)(**kwargs)

Create media_entry database table.

Contains:

entryId (PRIMARY KEY INTEGER),
[shotType](#) (TEXT NOT NULLABLE) an enum,
 time (TEXT NOT NULLABLE) format (yyyy-MM-ddTHH:mm:ss.zzz),
 illumination_type (TEXT NOT NULLABLE) an enum,
 iso (INTEGER NOT NULLABLE),
 apertureSize (REAL NOT NULLABLE),
 shutterSpeed (REAL NOT NULLABLE),
 whiteBalance (INTEGER NOT NULLABLE)

Method resolution order:

[MediaEntry](#)

sqlalchemy.orm.decl_api.Model

[flask_sqlalchemy.model.Model](#)[builtins.object](#)

Methods defined here:

__init__(self, **kwargs)

A simple constructor that allows initialization from kwargs.

Sets attributes on the constructed instance using the names and values in ``kwargs``.

Only keys that are present as attributes of the instance's class are allowed. These could be, for example, any mapped columns or relationships.

__repr__(self)

Return [MediaEntry](#).entryId when inserting into database.

Data descriptors defined here:

apertureSize**entryId****illuminationType****iso****shotType****shutterSpeed****time**

whiteBalance

Data and other attributes defined here:

```
__mapper__ = <Mapper at 0x7f2b6551e510; MediaEntry>
__table__ = Table('MediaEntry', MetaData(), Column('entryId'...table=<MediaEntry>, nullable=False), schema=None)
__tablename__ = 'MediaEntry'
```

Methods inherited from sqlalchemy.orm.decl_api.Model:

query

Data and other attributes inherited from sqlalchemy.orm.decl_api.Model:

```
__abstract__ = True
__fsa__ = <SQLAlchemy>
metadata = MetaData()
query_class = <class 'flask_sqlalchemy.query.Query'>
    SQLAlchemy :class:`~sqlalchemy.orm.query.Query` subclass with some extra methods
    useful for querying in a web application.

    This is the default query class for :attr:`~Model.query`.

    .. versionchanged:: 3.0
       Renamed to ``Query`` from ``BaseQuery``.
registry = <sqlalchemy.orm.decl_api.registry object>
```

Data descriptors inherited from [flask_sqlalchemy.model.Model](#):

```
__dict__
    dictionary for instance variables (if defined)
__weakref__
    list of weak references to the object (if defined)
```

Data and other attributes inherited from [flask_sqlalchemy.model.Model](#):

```
__annotations__ = {'__fsa__': 't.ClassVar[SQLAlchemy]', 'query': 't.ClassVar[Query]', 'query_class': 't.ClassVar[type[Query]]'}
```

class **MediaEntryInternalRepresentation**([builtins.object](#))

[MediaEntryInternalRepresentation](#)(entryId: int, [shotType](#): None, time: str, [illuminationType](#): None, iso: int, apertureSize: float, shutterSpeed: float, whiteBalance: i

Internal representation for [MediaEntry](#).

Methods defined here:

```
__init__(self, entryId: int, shotType: None, time: str, illuminationType: None, iso: int, apertureSize: float, shutterSpeed: float, whiteBalance: int)
    Construct MediaEntryInternalRepresentation.
```

Data descriptors defined here:

```
__dict__
    dictionary for instance variables (if defined)
__weakref__
    list of weak references to the object (if defined)
```

Data and other attributes defined here:

```
apertureSize = 0.0
entryId = 0
illuminationType = None
iso = 0
shotType = None
shutterSpeed = 0.0
```

time = "

whiteBalance = 0.0

```
class MediaMetadata(sqlalchemy.orm.decl_api.Model)
    MediaMetadata(**kwargs)
```

Create [MediaMetadata](#) database table.

Contains:

metadataId (PRIMARY KEY INTEGER),
 entryId (INTEGER, FOREIGN KEY ([MediaEntry.entryId](#))),
 leftCameraMedia (TEXT NOT NULLABLE),
 rightCameraMedia (TEXT NOT NULLABLE),
 time (TEXT NOT NULLABLE) format (yyyy-MM-ddTHH:mm:ss.zzzzzz),
 temperature (REAL NOT NULLABLE),
 pressure (REAL NOT NULLABLE),
 ph (REAL NOT NULLABLE),
 dissolvedOxygen (REAL NOT NULLABLE),

Method resolution order:

[MediaMetadata](#)
 sqlalchemy.orm.decl_api.Model
[flask_sqlalchemy.model.Model](#)
[builtins.object](#)

Methods defined here:

__init__(self, **kwargs)
 A simple constructor that allows initialization from kwargs.

Sets attributes on the constructed instance using the names and values in ``kwargs``.

Only keys that are present as attributes of the instance's class are allowed. These could be, for example, any mapped columns or relationships.

__repr__(self)
 Return [MediaMetadata.metadataId](#) when inserting into database.

Data descriptors defined here:

dissolvedOxygen

entryId

leftCameraMedia

metadataId

ph

pressure

rightCameraMedia

temperature

time

Data and other attributes defined here:

__mapper__ = <Mapper at 0x7f2b63cba7d0; MediaMetadata>

__table__ = Table('MediaMetadata', MetaData(), Column('metad...le=<MediaMetadata>, nullable=False), schema=None)

__tablename__ = 'MediaMetadata'

Methods inherited from sqlalchemy.orm.decl_api.Model:

query

Data and other attributes inherited from sqlalchemy.orm.decl_api.Model:

__abstract__ = True

`__fsa__` = <SQLAlchemy>

`metadata` = MetaData()

`query_class` = <class 'flask_sqlalchemy.query.Query'>
SQLAlchemy :class:`~sqlalchemy.orm.query.Query` subclass with some extra methods
useful for querying in a web application.

This is the default query class for :attr:`.Model.query`.

.. versionchanged:: 3.0
Renamed to ``Query`` from ``BaseQuery``.

`registry` = <sqlalchemy.orm.decl_api.registry object>

Data descriptors inherited from [flask_sqlalchemy.model.Model](#):

`__dict__`
dictionary for instance variables (if defined)

`__weakref__`
list of weak references to the object (if defined)

Data and other attributes inherited from [flask_sqlalchemy.model.Model](#):

`__annotations__` = {'__fsa__': 't.ClassVar[SQLAlchemy]', 'query': 't.ClassVar[Query]', 'query_class': 't.ClassVar[type[Query]]'}

class **MediaMetadataInternalRepresentation**([builtins.object](#))

[MediaMetadataInternalRepresentation](#)(metadataId: int, entryId: int, leftCameraMedia: str, rightCameraMedia: str, time: str, temperature: float, ph: float, dissolvedC

Internal representation for a [MediaMetadata](#) entry.

Methods defined here:

`__init__`(self, metadataId: int, entryId: int, leftCameraMedia: str, rightCameraMedia: str, time: str, temperature: float, ph: float, dissolvedOxygen: float)
Construct [MediaMetadataInternalRepresentation](#).

Data descriptors defined here:

`__dict__`
dictionary for instance variables (if defined)

`__weakref__`
list of weak references to the object (if defined)

Data and other attributes defined here:

`dissolvedOxygen` = 0.0

`entryId` = 0

`leftCameraMedia` = "

`metadataId` = 0

`ph` = 0.0

`pressure` = 0.0

`rightCameraMedia` = "

`temperature` = 0.0

`time` = "

class **illuminationType**([enum.Enum](#))

[illuminationType](#)(value, names=None, *, module=None, qualname=None, type=None, start=1, boundary=None)

Illumination Type.

Method resolution order:

[illuminationType](#)
[enum.Enum](#)
[builtins.object](#)

Data and other attributes defined here:

INFRARED = <illuminationType.INFRARED: 3>

NONE = <illuminationType.NONE: 1>

ULTRAVIOLET = <illuminationType.ULTRAVIOLET: 4>

VISIBLESPECTRUM = <illuminationType.VISIBLESPECTRUM: 2>

Data descriptors inherited from [enum.Enum](#):

name

The name of the Enum member.

value

The value of the Enum member.

Methods inherited from [enum.EnumType](#):

__contains__(member) from [enum.EnumType](#)

Return True if member is a member of this enum
raises TypeError if member is not an enum member

note: in 3.12 TypeError will no longer be raised, and True will also be
returned if member is the value of a member in this enum

__getitem__(name) from [enum.EnumType](#)

Return the member matching `name`.

__iter__() from [enum.EnumType](#)

Return members in definition order.

__len__() from [enum.EnumType](#)

Return the number of members (no aliases)

Readonly properties inherited from [enum.EnumType](#):

__members__

Returns a mapping of member name->value.

This mapping lists all enum members, including aliases. Note that this
is a read-only view of the internal mapping.

class **shotType**([enum.Enum](#))

[shotType](#)(value, names=None, *, module=None, qualname=None, type=None, start=1, boundary=None)

Shot Type for the type of study to be performed.

Method resolution order:

[shotType](#)
[enum.Enum](#)
[builtins.object](#)

Data and other attributes defined here:

BURST = <shotType.BURST: 2>

SINGLE = <shotType.SINGLE: 1>

TELESCOPIC = <shotType.TELESCOPIC: 3>

TIMELAPSE = <shotType.TIMELAPSE: 4>

VIDEO = <shotType.VIDEO: 5>

Data descriptors inherited from [enum.Enum](#):

name

The name of the Enum member.

value

The value of the Enum member.

Methods inherited from [enum.EnumType](#):

__contains__(member) from [enum.EnumType](#)

Return True if member is a member of this enum
raises TypeError if member is not an enum member

note: in 3.12 TypeError will no longer be raised, and True will also be returned if member is the value of a member in this enum

— **getitem__**(name) from [enum.EnumType](#)
Return the member matching `name`.

— **iter__**() from [enum.EnumType](#)
Return members in definition order.

— **len__**() from [enum.EnumType](#)
Return the number of members (no aliases)

Readonly properties inherited from [enum.EnumType](#):

— **members__**
Returns a mapping of member name->value.

This mapping lists all enum members, including aliases. Note that this is a read-only view of the internal mapping.

Functions

getAllMediaEntry(db) -> list[app.MediaEntryInternalRepresentation]
Get all Media Entries from a database.

getAllMediaEntryIDs(db)
Get all entryId from database.

getAllMediaMetadata(db) -> list[app.MediaMetadataInternalRepresentation]

getAllMediaMetadataId(db, entryId: int) -> list[app.MediaMetadataInternalRepresentation]

getCurrentTime() -> str
Return current time in 'yyyy-MM-ddTHH:mm:ss.zzz' format.

index()
Return index.html to the / route.

insertMediaEntry(db, shotType, illuminationType, iso, apertureSize, shutterSpeed, whiteBalance)
Insert a [MediaEntry](#) into the database.

insertMediaMetadata(db, entryId, leftCameraMedia, rightCameraMedia, temperature, pressure, ph, dissolvedOxygen)
Insert [MediaMetadata](#) entry into database.

Data

CONTINUOUS = <EnumCheck.CONTINUOUS: 'no skipped integer values'>

UNIQUE = <EnumCheck.UNIQUE: 'one name per value'>

app = <Flask 'app'>

basedir = '/home/uwu/Github_Repos/University/M.O.S.I.S/M.O.S.I.S_Host_Software/App'

db = <SQLAlchemy>